

HARDWARE REQUIREMENTS

This chapter defines the system requirements for ENDOCLICK / +.

QlikView	32-bit (x86)	64-bit (x64)
Operating System	Windows 2000 Pro SP4 ¹ Windows 2000 Server SP4 ¹ Windows XP Windows Vista Windows Server 2003 Windows Server 2008	Windows XP Professional x64 Edition Windows Vista x64 Windows Server 2003 x64 Edition Windows Server 2008 x64 Edition
Processor	Pentium 4 compatible or higher recommended	x64 compatible processors
Memory	Depending on size of applications; 512 MB or more recommended	Depending on size of applications; 2 GB or higher recommended
Disk space	150 MB total required to install	150 MB total required to install
Security	Windows Integrated Security or custom authentication ²	Windows Integrated Security or custom authentication ²

QlikView Web Clients	QlikView Analyzer for Microsoft Internet Explorer	QlikView Analyzer Java Client	QlikView Analyzer AJAX Client (Zero Footprint)
Browser	Internet Explorer 6.0 or 7.0 (Windows)	Internet Explorer 6.0 or 7.0 (Windows) Firefox 2.0 (Windows, Mac, Unix) Safari 2.0 or above (Windows, Mac)	Internet Explorer 6.0 or 7.0 (Windows) Firefox 2.0 (Windows, Mac, Unix) Safari 2.0 or above (Windows, Mac)
Java Virtual Machine	Not applicable	Sun Java Runtime 5.0 or above	Not applicable
QlikView Server	Small Business and Enterprise Edition	Small Business and Enterprise Edition	Enterprise Edition
Security	Windows Integrated Security or custom authentication ²	Windows Integrated Security or custom authentication ²	Windows Integrated Security or custom authentication ²

QlikView Server	32-bit (x86)	64-bit (x64)
Platform	Windows XP Professional ³ Windows Vista ³ Windows Server 2003 Windows Server 2008	Windows XP Professional x64 Edition ³ Windows Vista x64 ³ Windows Server 2003 x64 Edition Windows Server 2008 x64 Edition
Processor	Pentium 4 compatible or higher recommended	x64 compatible processors
Memory	Depending on size of applications; 1 GB or more recommended	Depending on size of applications; 4 GB or more recommended
Disk space	450 MB total required to install	450 MB total required to install
Security	Windows Integrated Security or custom authentication	Windows Integrated Security or custom authentication
Web server	QlikView http Server or Microsoft IIS 6.0 or 7.0 ⁴	QlikView http Server or Microsoft IIS 6.0 or 7.0 ⁴

QlikView Publisher	32-bit (x86)	64-bit (x64)
Platform	Windows XP Professional ³ Windows Vista ³ Windows Server 2003 Windows Server 2008	Windows XP Professional x64 Edition ³ Windows Vista x64 ³ Windows Server 2003 x64 Edition Windows Server 2008 x64 Edition
Processor	Pentium 4 compatible or higher recommended	x64 compatible processors
Memory	Depending on size of applications; 512 MB or more recommended	Depending on size of applications; 2 GB or more recommended
Disk space	250 MB total required to install	250 MB total required to install
Database	SQL Server 2000 ⁵ , SQL Server 2005 ⁵ , Native XML	SQL Server 2000 ⁵ , SQL Server 2005 ⁵ , Native XML
Web server	QlikView http Server or Microsoft IIS 6.0 or 7.0	QlikView http Server or Microsoft IIS 6.0 or 7.0
.NET	2.0.50727	2.0.50727
Security	Windows Integrated Security or custom authentication ⁶	Windows Integrated Security or custom authentication ⁶
Browser	Internet Explorer 6.0 or 7.0	Internet Explorer 6.0 or 7.0