

OLYMPUS

Your Vision, Our Future

OLYMPUS STONE MANAGEMENT EZDILATE BALLOON DILATORS

Precise Access and Dilation for PCNL and URS.



PRECISE VISUALIZATION FOR CONTROLLED PERCUTANEOUS AND URETEROSCOPIC ACCESS

EZDilate balloon dilators are available in multiple sizes to allow for safe ureter dilation for stricture relief as well as safe access during percutaneous surgery with minimized patient trauma.


- Sheath is tapered at the distal end for **smooth, atraumatic insertion** into the patient's anatomy even when there is challenging surrounding tissue.
- **Safe and precise visualization** thanks to fully radiopaque tip ensuring access to the desired location prior to inflation.
- **Controlled access** thanks to accurate and predictable **single-step dilation**. Specialized hydrophilic and silicone coatings on the balloon promote smooth, low-friction positioning.

Nephro-EZDilate Balloon Dilator

Precise visualization due to radiopaque tip

Controlled single-step dilation


Smooth and atraumatic insertion



10787

Uro-EZDilate Balloon Dilator

Designed to ensure stricture relief and gentle ureteral dilation, the balloon is available in 4 and 10 cm and from 12 to 18 Fr to accommodate any procedural requirements.



10973

Balloon Inflation Device

The high-pressure, disposable inflation device is included with all EZDilate balloon dilator kits and is also available separately. Its intuitive design includes a simple plunger-engagement switch that allows for easy, precise inflation as well as rapid deflation.



10795

Folding Tool

Integrated folding tool ensures small profile upon reinsertion.



10783

EZDILATE BALLOON DILATORS

Ordering Information

Nephro-EZDilate Balloon Dilator – Burst Pressure 17 atm

Kit: Balloon Dilator with 25 cc Inflation Device

Catalogue No.	Description	Dimensions
EGBPCN0815K	PCNL balloon and inflation device kit	8 mm (24 Fr) × 15 cm
EGBPCN1015K	PCNL balloon and inflation device kit	10 mm (30 Fr) × 15 cm

Individual Balloon Dilators

Catalogue No.	Description	Dimensions
EGBPCN0815	PCNL balloon	8 mm (24 Fr) × 15 cm
EGBPCN1015	PCNL balloon	10 mm (30 Fr) × 15 cm
EGBINF2530	Balloon inflation device	-

Uro-EZDilate Balloon Dilator – Burst Pressure 20 atm

Kit: Balloon Dilator with 25 cc Inflation Device

Catalogue No.	Description	Dimensions
EGBURS0404K	URS balloon and inflation device kit	4 mm (12 Fr) × 4 cm
EGBURS0410K	URS balloon and inflation device kit	4 mm (12 Fr) × 10 cm
EGBURS0504K	URS balloon and inflation device kit	5 mm (15 Fr) × 4 cm
EGBURS0510K	URS balloon and inflation device kit	5 mm (15 Fr) × 10 cm
EGBURS0604K	URS balloon and inflation device kit	6 mm (18 Fr) × 4 cm
EGBURS0610K	URS balloon and inflation device kit	6 mm (18 Fr) × 10 cm

Individual Balloon Dilators

Catalogue No.	Description	Dimensions
EGBURS0404	URS balloon	4 mm (12 Fr) × 4 cm
EGBURS0410	URS balloon	4 mm (12 Fr) × 10 cm
EGBURS0504	URS balloon	5 mm (15 Fr) × 4 cm
EGBURS0510	URS balloon	5 mm (15 Fr) × 10 cm
EGBURS0604	URS balloon	6 mm (18 Fr) × 4 cm
EGBURS0610	URS balloon	6 mm (18 Fr) × 10 cm
EGBINF2530	Balloon inflation device	-

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG
Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14–18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com