

## TC-C3 Compact Trolley



# Compact Trolley

## TC-C3

The TC-C3 Compact Trolley is designed as a compact trolley to be used in conjunction with Optera CV-170.

### Compact

The compact design means the TC-C3 is ideal for use in the outpatient or doctor's office where space is limited.

### Easy to Move

The Compact Trolley is easy to move and position due to its light weight and ergonomically designed handles at the front of the top tray.

### Easy to Clean

The TC-C3 is easy to clean due to limited dirt traps and no exposed fittings.

### Monitor Support

The TC-C3 is offered with a monitor mount that supports monitors up to 26" with unrestricted swivel and tilt movement, allowing good viewing positions of the procedure image.

### Large Top Tray

The top tray provides a clear working space at a comfortable height. The increased range of accessories enables the trolley to be utilized for a range of procedures.

### Integrated Separation Transformer

The TC-C3 can be supplied with a separating transformer that is placed under the base, this means increased shelf space for equipment loading.



# Accessories

## TC-C3

### MAJ-2146

Sliding Keyboard Tray



### MAJ-1684

Scope Pole Kit



### MAJ-2148

Rear Panel



### MAJ-1686

Drawer Unit



### MAJ-1685

IV Bag Hanger



### MAJ-1683

LCD Monitor Mount



### MAJ-1656

Footswitch Holder



### MAJ-1639

CO<sub>2</sub> Cylinder Holder (140 mm)



### MAJ-1665

Camera Head Holder



### MAJ-2145

ENT Scope Hanger



### MAJ-1658

Universal Storage Container



### MAJ-110

Jar Mounting Kit TC



# TC-C3

## Standard sets include:

TC-C3 Compact Trolley plus items marked (•)

Accessories available as an optional extra (o)

			Region	Set 1	Set 2	Set 3	Set 4
MAJ	Part	Part	EUROPE				
Description	Description	Number	UK				
			ROW				
MAJ-2146	Sliding Keyboard Tray	K10027573	o	•	•	•	•
MAJ-1684	Scope Pole Kit	K10027281	o	•	•		
MAJ-2148	Rear Panel	K10027574	o				
	Transformer			•		•	
MAJ-1686	Drawer Unit	K10027569	o				
MAJ-1685	IV Bag Hanger	K10027568	o				
MAJ-1683	LCD Monitor Mount	K10026888	o	•	•	•	
MAJ-1656	Footswitch Holder	K10021790	o				
MAJ-1639	CO <sub>2</sub> Cylinder Holder (140 mm)	K10021041	o				
MAJ-1665	Camera Head Holder	K10028704	o				
MAJ-2145	ENT Scope Hanger	K10026497	o				
MAJ-1658	Universal Storage Container	K10021792	o				
MAJ-110	Jar Mounting Kit TC	K7503556	o				
MAJ-1654	Equipotential Terminal Strip	K10021352	o				
	Blister Pack Toolkit		o	•	•	•	•
	IEC Lead Set		o	•		•	

## Specifications

<b>Dimensions</b>	1,110 mm (H) × 563.9 mm (D) × 489.5 mm (W)
	Height to top tray work surface – 1,100 mm
	Height with 26" LCD monitor in highest position – 1,945 mm
<b>Weight</b>	Maximum 82 kg (unladen) including fitted transformer

## Voltage

The TC-C3 can be supplied with a separating transformer based on one of the following voltage options:

MAJ-2205	Transformer 110-120 V
MAJ-2206	Transformer 220-240 V

## Trolley Construction

<b>Columns</b>	Extruded profiled aluminum
<b>Shelves</b>	Mild steel
<b>Top Tray</b>	Injection-moulded ABS

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

**OLYMPUS**

**OLYMPUS EUROPA SE & CO. KG**  
 Postbox 10 49 08, 20034 Hamburg, Germany  
 Wendenstrasse 14-18, 20097 Hamburg, Germany  
 Phone: +49 40 23773-0, Fax: +49 40 233765  
 www.olympus-europa.com