

OLYMPUS

Your Vision, Our Future

IMAGING EQUIPMENT

EndoLED

Miniature Light Source

**MOBILITY
WITHOUT
COMPROMISE**



Original Size

EndoLED LIGHT SOURCE


FLEXIBILITY WITHOUT COMPROMISE IN BRIGHTNESS AND CONVENIENCE

The Olympus EndoLED combines the latest LED technology delivering unprecedented white light with an excellent colour reproduction.


Each type of the high power EndoLED uses a product specific light focusing technique to achieve outstanding brightness.

High quality material minimizes reprocessing concerns by allowing even Sterrad sterilization including Olympus recommended batteries.


- Colour temperature like Xenon light thanks to careful LED selection
- High brightness over the entire LED lifetime of 100,000 hrs
- Convenient battery concept using standard AA rechargeable batteries and a battery status indicator
- Ergonomic design for rigid and flexible endoscopes for more comfort
- Highest reprocessing standards: fully immersible, washer and disinfectant compatible, sterilizable with Sterrad



e.g HD Laryngoscope in ENT



e.g CYF-5 in Urology



BATTERY INDICATOR

For more safety, the battery status indicator shows green for fully charged, yellow for partially discharged or flashing yellow for fully discharged.

EndoLED LIGHT SOURCE

Technical Data


WA91500A Light source "EndoLED",
for rigid endoscopes

for rechargeable batteries,
with battery status indicator
(last up to 3 hrs),
includes 2 adapters O0332




WA91502A Light source "EndoLED",
for flexible endoscopes


for rechargeable batteries,
with battery status indicator
(last up to 3 hrs)



WA91505A Battery,
AA type,
rechargeable,
4 pcs.



System Chart



OLYMPUS

Published by
OLYMPUS WINTER & IBE GMBH

Kuehnstraße 61, 22045 Hamburg, Germany

Distributed by
OLYMPUS MEDICAL SYSTEMS CORPORATION
OLYMPUS MEDICAL SYSTEMS EUROPA GMBH
OLYMPUS AMERICA INC.
KEYMED (MEDICAL & INDUSTRIAL EQUIPMENT) LTD.
OLYMPUS SINGAPORE PTE LTD.
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
OLYMPUS AUSTRALIA PTY. LTD.
OLYMPUS LATIN AMERICA INC.
OLYMPUS KOREA CO., LTD.
OLYMPUS (BEIJING) SALES & SERVICE, CO., LTD.

2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
Wendenstraße 14-18, 20097 Hamburg, Germany
3500 Corporate Parkway, P.O. Box 610, Center Valley, PA 18034-0610, U.S.A.
KeyMed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, United Kingdom
491B, River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373, Singapore
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
31 Gilby Road, Mount Waverley, Victoria 3149, Australia
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126-2097, U.S.A.
4F, Gyeongam Bldg., 157-27, Samseong-Dong, Kangnam-Gu, Seoul 135-090, Korea
R1202, NCI Tower, A12 Jianguomenwai Dajie, Chaoyangqu Beijing 100022, China